

ChronoGate

Standard Operating Procedure — interactive time-gating, lifetime & phasor FLIM analysis

Version 0.18.0

2026-07-21

PicoQuant .ptu · Becker & Hickl .sdt

1. Purpose & scope
2. What ChronoGate does
3. Install & launch
4. Opening data
5. The workspace
6. Setting the gate and t_0
7. Display controls
8. Spatial binning
9. Lifetime & phasor analysis
 - 9.1 Two-gate RLD
 - 9.2 Phasor
10. IRF reconvolution fit (rigorous lifetime)
11. Selecting pixels & regions
12. Exporting & the one-page report
13. Troubleshooting
14. Revision history

1. Purpose & scope

This procedure describes how to operate **ChronoGate**, an interactive time-gating, rapid-lifetime and phasor viewer for PicoQuant `.ptu` (and Becker & Hickl `.sdt`) FLIM data. It covers loading data, setting the time gate, producing lifetime and phasor maps, running a rigorous IRF-reconvolution fit, and exporting results.

It is a practical operating guide, not a validation report. Absolute-lifetime work should use a measured IRF (see §10) and be cross-checked against a reference of known lifetime.

2. What ChronoGate does

A pulsed laser excites the sample; for every detected photon the hardware records the delay since the last pulse (the **microtime**). The per-pixel histogram of those delays is the **fluorescence decay**. ChronoGate reconstructs that decay cube and lets you:

- **Time-gate** — keep photons in a chosen microtime window and sum them into an image (fit-free lifetime *contrast*).
- **RLD** — read an apparent lifetime from two gates (fit-free, per pixel).
- **Phasor** — plot each pixel's decay as a point on the universal semicircle.
- **IRF fit** — recover a rigorous, IRF-deconvolved lifetime by reconvolution.




Figure 1. The end-to-end operating flow: load → gate → choose an analysis → export.

3. Install & launch

ChronoGate runs from source in a Python **3.12+** environment.

macOS (double-click launcher). Keep `ChronoGate.app` and `ChronoGate.command` together in the project folder and double-click `ChronoGate.app`. The first launch builds a virtual environment from `pyproject.toml` (a one-time download); later launches are instant. Changing dependencies triggers an automatic rebuild.

From a terminal (any OS).

```
python -m chronogate          # opens a file picker
python -m chronogate path/to/file.ptu
python -m chronogate path/to/file.sdt
```

If opening a `.sdt` reports that `sdtfile` is missing, or an IRF fit reports `scipy` is missing, your environment predates those dependencies — relaunch so the launcher rebuilds it, or run `pip install -e .`

4. Opening data

- **File > Open (Ctrl+O)** — choose a `.ptu` or `.sdt` file.
- **File > Open folder (stack)** — load a numbered z-series (e.g. `FLIM_stack_z1.ptu ... z65.ptu`); step through planes with the **z-slice** slider in the File panel.
- **Channel** — for multi-detector files, pick the detector channel.

5. The workspace




Figure 2. The main window. **1** File/layer (open, z-slice, channel, export, settings). **2** Gate (start, end, t0). **3** Display (threshold, noise floor, colormap, scale lock). **4** Binning. **5** Lifetime (RLD) controls. **6** Per-pixel decay / selection tools. The left plot is the summed decay with the gate shaded; the right plot is the current image.

Switch the image with the mode buttons in the toolbar: **Intensity** (I), **Lifetime (RLD)** (T), **Phasor** (P).

6. Setting the gate and t0



Figure 3. Gate panel: **1** gate start, **2** gate end (both in ns; you can also drag the shaded span on the decay plot), **3** t0 (pulse reference), **4** auto (set t0 to the smoothed decay peak), **5** reset (set t0 to its default, 0 ns).

1. Drag the shaded band on the decay plot, or type **start/end** in ns.
2. Set **t0** to the pulse arrival: click **auto** for the detected peak, or type a value. The **reset** (↺) button returns t0 to 0 ns.
3. The image updates live as you move the gate.

7. Display controls




Figure 4. Display panel: **1** min photons/px (blank dim pixels), **2** noise floor (counts, read on the summed decay), **3** floor auto (robust baseline), **4** floor reset (0 — no subtraction), **5** colormap.

- **min photons/px** masks pixels too dim to trust.
- **noise floor** subtracts a flat background so a pedestal doesn't bias the lifetime. **auto** picks a robust baseline; **reset** turns subtraction off.
- **lock scale** freezes the colour range so planes/frames are comparable.

8. Spatial binning



Figure 5. Binning panel: 1 bin factor, 2 target photons/px, 3 *Auto* (suggest a factor from photon statistics), 4 *reset* (1×, off).

FLIM pixels are often photon-starved. Binning pools each pixel's neighbourhood ($b \times b$) for cleaner decays and lifetimes. Click **Auto** to reach the target photon budget, or set the factor by hand. **reset** (↺) returns to 1× (off).

Binning is the fastest way to make a per-pixel lifetime *map* usable on sparse data — see §9 and §10.

9. Lifetime & phasor analysis

9.1 Two-gate RLD

Press **T** for Lifetime mode. RLD reads an apparent lifetime from two equal-width gates, per pixel, with no fitting:




Figure 6. Two-gate RLD. With equal-width gates separated by Δt , the width/amplitude factors cancel, giving $\tau = \Delta t / \ln(N_A/N_B)$. The shot-noise uncertainty σ_τ is reported for a pooled region.

1. In the Lifetime panel choose which gate you are editing (**A** early / **B** late) and position the two gates on the decay.
2. Keep the gates **equal width** (the panel warns if not) — RLD assumes it.
3. The τ image appears on the right. Select a region (§11) and open the one-page report (§12) to get the pooled $\tau \pm \sigma$ (shot-noise uncertainty).

9.2 Phasor

Press **P**. Each pixel's decay becomes a point (g, s); a single lifetime lands on the universal semicircle, mixtures fall inside it. Use **Calibrate phasor from reference...** (View menu) with a dye of known lifetime to remove the instrument/t0 offset, and the **2nd harmonic** toggle to spread short lifetimes apart.

10. IRF reconvolution fit (rigorous lifetime)

For an IRF-deconvolved lifetime (rather than fit-free contrast), use **View** ▸ **IRF lifetime fit...**

The screenshot shows the 'IRF fit dialog' with the following settings and callouts:

- Instrument response (IRF)**:
 - 1. **Gaussian model** (selected)
 - centre: 0.500 ns
 - FWHM: 0.291 ns
 - 2. **Measured IRF file...** (selected)
 - file: no IRF file selected
 - Browse...
- Model**:
 - 3. **components**: mono (1 exponential)
 - 4. **objective**: Poisson MLE (recommended)
 - 5. **photon threshold**: 50
- Fit scope**:
 - 6. **Selected region — one fit (fast)** (selected)
 - Whole image — per-pixel τ -map (slow)

Reconvolution τ is IRF-deconvolved — distinct from the fit-free RLD τ and phasor. The report labels each by method.

OK Cancel

Figure 7. IRF fit dialog: 1 Gaussian IRF model, 2 measured IRF file, 3 components (mono/bi), 4 objective (Poisson MLE / weighted χ^2), 5 photon threshold, 6 whole-image τ -map (vs. a region fit).

1. **Choose the IRF.** For a first look, a **Gaussian** model is fine. For **absolute** lifetimes, load a **measured IRF** file (a scatter/reflection or reference-dye acquisition) — reconvolution is only as good as its IRF.
2. **Model** — mono for one lifetime, bi for two.
3. **Objective** — Poisson MLE (recommended for low counts) or weighted χ^2 .

4. **Scope** — a **region fit** (select a bright ROI first) gives one trustworthy $\tau \pm \sigma$; a **whole-image τ -map** fits every pixel above the photon threshold.
5. Read the result: **reduced $\chi^2 \approx 1$** means a good fit. A large χ^2 (e.g. tens or hundreds) means the model is wrong — try a measured IRF, or bi-exponential (see Troubleshooting). A component reported as **$\sigma = n/a$ (unidentifiable)** is an honest flag that the data cannot separate it.

For a per-pixel map on sparse data, **bin first** (§8) so pixels clear the photon threshold, then run the map. The default threshold adapts to the image.

11. Selecting pixels & regions

- **Drag a box** on the image for a rectangular ROI; **lasso** on the phasor plot to select a population; or use the **Pixel list** (View ► Pixel list) to pick by metric.
- Hover to preview a pixel's decay; click to lock it; **Pin** to keep several for comparison.
- A selection travels with every export and scopes the report's statistics.

12. Exporting & the one-page report




Figure 8. Export dialog: 1 raster TIFF (raw values, for ImageJ), 2 one-page report (PNG+PDF with the headline number, decay, image and provenance), 3 restrict to the current selection, 4 *Export & open in Fiji*.

- **File ► Export** (**Ctrl+E**) — choose artefacts (raw TIFF, colormapped PNG, decay CSV, per-pixel table, selection files) and a folder. A content-hashed **provenance JSON** is always written.
- **File ► Export report** (**Ctrl+R**) — the one-page summary; in Lifetime mode it includes the pooled $\tau \pm \sigma$ for the selection.
- **Restrict to selection** masks everything outside the selected pixels.

- **Export & open in Fiji** writes a raw TIFF plus an ImageJ macro (range, LUT, ROI) and launches Fiji on it (set the Fiji path in Preferences first).

13. Troubleshooting

Symptom	Cause	Fix
τ -map fits 0 px	photon threshold above every pixel's count	bin more (§8), lower the threshold, or do a region fit
IRF fit reduced $\chi^2 \approx 100s$	wrong/guessed IRF, or model too simple	load a measured IRF; try bi-exponential; check the Gaussian centre/FWHM
σ shows n/a (unidentifiable)	near-degenerate fit (e.g. two equal τ)	expected/honest — simplify the model or accept the point estimate
<code>.sdt</code> won't open (<code>sdtfile</code> missing)	environment predates the dependency	relaunch (auto-rebuild) or <code>pip install -e .</code>
RLD τ looks wrong	unequal gates, or no background subtraction	make gates equal width; set the noise floor (§7)
app "unidentified developer" warning	installers are unsigned	right-click ▶ Open (macOS) / More info ▶ Run (Windows)

14. Revision history

The SOP is regenerated from the app: edit `docs/sop/sop.md` and `docs/sop/figures.py`, then run `python docs/sop/build.py`.

Version	Date	Notes
0.18.0	2026-07-21	Generated from ChronoGate 0.18.0.